## **Bouffant Cap**

## Features

Made of Spunbonded polypropylene, and sewn with elastic band, with its class one flammability rating and latex free features, it gives your head great protection and comfort. There are size 18",21" and 24" available, which fit most people. It's widely used for cleaning, food service, healthcare, laboratory and homecare.

- Latex free
- Suitable for isolation and basic protection from bacteria and particulate
- Specialized type or design for particular occasions
- Practical barrier for wide range of applications
- Soft and light weight
- Good fit, feel and performance



## Applications



Medical Purpose / Examination



Industrial purpose / PPE



Laboratory



Healthcare and nursing



General housekeeping



IT Industry











SKU	Size	Color	Material	Package	Carton Dimension
BC101002	21"	White	SPP 10GSM	100pcs/pk,10pks/ctn	340*270*270mm
BC101001	21"	Blue	SPP 10GSM	100pcs/pk,10pks/ctn	340*270*270mm
BC141002	21"	White	SPP 14GSM	100pcs/pk,10pks/ctn	340*270*270mm
BC141001	21"	Blue	SPP 14GSM	100pcs/pk,10pks/ctn	340*270*270mm

